




DLF

**DYDAKTYCZNE
LABORATORIUM
FIZYCZNE**

Instytut Fizyki Doświadczalnej
Wydział Matematyki, Fizyki i Informatyki
UNIwersytet Gdański


Experiment 26

Studying ferromagnetic properties based on hysteresis loops



I. Background theory.

1. Magnetisation, magnetic permeability, magnetic susceptibility.
2. Types of magnetism.
3. Ferromagnetism:
 - a) ferromagnetic configurations;
 - b) Curie – Weiss law;
 - c) energy: exchange, anisotropy, magnetic;
 - d) domain structure in solids;
 - e) primary magnetization curve of ferromagnetic materials;
 - f) hysteresis loop.
4. Other types of magnetic moments.



Picture 1. Experimental setup for studying ferromagnetic properties: 1 – oscilloscope; 2 – power dial; 3 – phase shifter; 4 – ferromagnetic coil; 5 – universal multimeter; 6 – thermostat; 7 – transformer.

II. Experimental tasks.

1. Familiarise yourself with the setup shown in *Picture 1* and with the component connection diagram in *Figures 2 and 3* and in the *Appendix*.
2. Determine the initial magnetization curve of iron at room temperature.
To do this, adjust the voltage U across the coils from 6 to 24 V.
Turn on the transformer (7 in *Picture 1*) and select regulated voltage V_{\sim} . Adjust the transformer settings as indicated by the multimeter (5 in *Picture 1*).
Read off the corresponding hysteresis loop parameters from the oscilloscope screen (1 in *Picture 1*) for each coil voltage.
3. Determine the dependence of the area under the hysteresis loop P on the magnetic field strength H for each voltage U as in step 2 above.

- Determine the dependence of magnetization M of the iron test sample on temperature T .
Measure the size of the hysteresis loop at constant magnetic field strength for temperatures ranging from room temperature to 150 °C.
To do this, set the coil voltage on the meter (5 in *Picture 1*) to about 10 V.
Measure the height of the hysteresis loop at room temperature.
Turn on the power controller and thermostat (2 and 6 in *Picture 1*). Set the heating knob on the power controller to position 10.
After reaching 150 °C, turn off the controller and start measuring the hysteresis loop (while the sample cools).
- Present the results of steps 2 – 4 graphically. Make the following plots:
 $M = f(H)$, $P = f(H)$, $1/M = f(T)$.
Interpret the curves obtained.
- Estimate the Curie – Weiss temperature for the iron sample.

III. Apparatus.

- Oscilloscope.
- Power controller.
- Phase shifter.
- Set of coils with iron sample.
- Universal multimeter.
- Thermostat.
- Transformer.

IV. Literature.

- R.P. Feynman, R. Leighton, M.Sands – *“The Feynman Lectures on Physics”*, Wesley 2005.
- Ch. Kittel – *“Introduction to Solid State Physics”*, Wiley, 2004.
- H. A. Enge, M.R. Wehr, J.A. Richards – *“Introduction to Atomic Physics”*, Wesley, 1981.

Appendix

Experimental setup wiring diagrams




Figure 2. General wiring diagram of the experimental setup.




Figure 3. Wiring diagram of the set of coils with phase shifter.